

Effective Pressure Area Care

BARIATRIC

EHOB®'s dedicated range of Bariatric Pressure Area Care products


Ventilation Holes:

- aiding the even distribution of pressure across the patient's body surface area.
- provides effective ventilation and optimises micro-climate.
- minimises risk of skin maceration as body fluids are drawn away from the patient's skin.

Reinforced handwells



4250A Bariatric Cushion
(LAD pump sold separately)



Indications for use:

- to aid in the prevention of pressure ulcers
- can be used for all levels of risk
- to aid in the treatment of all categories of pressure ulcers (I - IV)



4010A Bariatric Overlay Mattress
including LAD pump

Material:

High grade, antimicrobial, TPU material with optimal envelopment and immersion properties to ensure clinical efficacy and patient comfort.

EHOB® products have been clinically evaluated as part of effective protocol for the prevention and treatment of pressure ulcers.



EHOB® has provided effective solutions for prevention and treatment of pressure ulcers for more than 30 years.

We offer an extensive range of highly effective pressure relieving products to fit your specific needs.

EHOB® Bariatric Heel Protectors



Static Air

4420A Heel elevator TPU Custom XL

- Complete heel elevation
- Fully adjustable
- Suitable for multipatient use
- Silicone Buttons
- Optional Extender Straps



Fibre Filled

TRUVUEXW060 TruVue Heel Protector XL

- Complete heel elevation
- Fully adjustable
- Machine washable
- Special mix of fibres to ensure optimal pressure redistribution
- Optional Antirotation Wedge

EHOB® Bariatric Product Range:

Order number	Description	Size (inflated)	Weight Limit
4010A	EHOB Bariatric Overlay Mattress	192 cm x 120 cm x 8 cm	360 kg
4250A	EHOB Bariatric Cushion	71 cm x 56 cm x 5 cm	320 kg
4420A	EHOB Heel Elevator TPU Custom XL	Calf Size 46 - 64 cm	N.A.
TRUVUEXW060	TruVue Heel Protector - Extra Large	Calf Size 46 - 64 cm	N.A.

520 LAD hand pump sold separately

Can be reused following local policy guidelines for cleaning and disinfection. All EHOB products in the Static Air range incorporate an antimicrobial agent that acts as an effective barrier against the following microorganisms: Methicillin-Resistant Staphylococcus Aureus (MRSA); Clostridium difficile; Escherichia coli; Klebsiella pneumonia; Pseudomonas aeruginosa; Salmonella choleraesuis; Candida albicans; Enterococcus faecium; Aspergillus niger and Bacillus subtilis.



H5 Healthcare
2 /15 Seasands Drive Redhead NSW 2290
Tel.: 1300 727 631 - info@h5healthcare.com.au
www.h5healthcare.com.au