




H5 Healthcare

Effective Pressure Area Care Heel Cushion

High grade, antimicrobial TPU material with optimal *envelopment and immersion* properties to ensure clinical efficacy and patient comfort



Static Air Design

- Promotes even distribution of weight to prevent and treat pressure injuries
- provides effective ventilation
- comfortably supports by flotation, minimising tissue distortion by offloading the bony structures of the feet

- ✓ Easy to use
- ✓ Latex-free
- ✓ 6 month warranty

Code 4310 86cm x 43cm x 9cm

10.5 strokes to inflate

The secure strap can be trimmed to your required length

The Heel Cushion should be placed at the foot end of the bed using the secure strap. The inflation valve is to be placed face down against the mattress so that it does not irritate the patients feet.

Can be reused following local policy guidelines for cleaning and disinfection.

All EHOB products in the Static Air range incorporate an antimicrobial agent

that acts as an effective barrier against the following microorganisms: Methicillin-Resistant Staphylococcus Aureus (MRSA); Clostridium difficile; Escherichia coli; Klebsiella pneumonia; Pseudomonas aeruginose; Salmonella cholerausuis; Candida albicans; Enterococcus faecium; Aspergillus niger and Bacillus subtilis;



H5 Healthcare

H5 Healthcare
2/15 Seasands Drive Redhead NSW 2290
Tel: 1300 727 631 sales@h5healthcare.com.au
www.h5healthcare.com.au

EHOB[®]
prevention • treatment • solutions