



WAFFLE® Product Range


The WAFFLE® Products are made from high grade antimicrobial TPU material with optimal envelopment and immersion properties to ensure clinical efficacy and patient comfort.

VENTILATION HOLES:

- Aiding the even distribution of pressure across the patients body surface area.
- Provides effective ventilation and optimises micro-climate
- Minimises risk of skin maceration as body fluids are drawn away from the patient skin.


INDICATIONS FOR USE:

- To aid in the prevention of pressure ulcers
- Can be used for all levels of risk
- To aid in the treatment of all categories of pressure ulcers (I-IV)



BENEFITS

- » Highly cost effective
- » Easy to setup
- » Latex Free
- » Up to 372kg weight limit depending on product
- » Antimicrobial agent incorporated which is effective against a wide range of microorganisms including MRSA, C.Diff and E.Coli
- » Low profile
- » 6 month warranty



WAFFLE® OVERLAY

- The 'dog bone' handwells facilitate easy log rolling, mobilising and patient transfers
- Corner Straps secure and centre the overlay on the bed
- Easy Inflation
- Clear valve for inflation and adjustments
- Red CPR valve for rapid deflation

WAFFLE® HEEL PROTECTOR

- The Static Air® design ensures offloading of the heel and effective protection of the foot, ankle and lower limb.
- Easy access to the heel for dressing changes and visual inspection of the heel and foot plantar surface. Heel is continuously offloaded.



TruVue® Heel Protector

The TruVue® Heel Protector is filled with a special mix of fibres which ensures an optimal micro-climate. The Achilles Cushion lifts the heel off the support surface.

- The TruVue® Heel Protector is built up around the foot and lower limb for optimal levels of comfort and support.
- The opening sole allows easy inspection of the heel and plantar surface of the foot. The extra long sole provides protection of the toes
- Adjustable straps to keep the foot in an optimal position and minimise the risk of foot drop and the development of ankle contractures
- Protection against shear forces
- Internal/external rotation minimised with use of the wedge



BENEFITS

- » Rigid adjustable anti foot-drop straps
- » Anti-shearing barrier pad
- » Anti-rotation wedge (optional)

CODE	DESCRIPTION	SIZE	WEIGHT LIMIT
210WC	WAFFLE® Cushion Original	43x43 cm	136kg
228WC	WAFFLE® Cushion Standard	48x48 cm	136kg
4250	WAFFLE® Cushion Bariatric	71x56 cm	320kg
4310	WAFFLE® Heel Immersion Cushion	86x43x9 cm	136kg
4000	WAFFLE® Overlay	193x86x8 cm	270kg
4010	WAFFLE® Overlay Bariatric	193x119x8 cm	372kg
4290	WAFFLE® Trolley Overlay	178x51x7 cm	230kg
4440	WAFFLE® Heel Protector	One Size – Velcro Straps	
540MAD	Large Hand Pump	Suitable 4000 4010 4290 4310	
520LAD	Small Hand Pump	Suitable 4230 4240 4250 4440	
TRUVUEPW	TruVue® Petite with Wedge	Petite	
TRUVUESW	TruVue® Standard with Wedge	Standard	
TRUVUEXW	TruVue® Extra Large with Wedge	Extra Large	
TRUVUEPT	TruVue® Petite	Petite	
TRUVUESD	TruVue® Standard	Standard	
TRUVUEXL	TruVue® Extra Large	Extra Large	

